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- [1] Hsu W, Lesniak MS, Tyler B, Brem H. Local delivery of interleukin-2 and adriamycin is synergistic in the treatment of experimental malignant glioma. *J Neurooncol* 2005;74:135–140.
- [2] Lesniak MS, Gabikian P, Tyler BM, Pardoll DM, Brem H. Dexamethasone mediated inhibition of local IL-2 immuno-therapy is dose dependent in experimental brain tumors. *J Neurooncol* 2004;70:23–28.
- [3] Lesniak MS, Upadhyay U, Goodwin R, Tyler B, Brem H. Local delivery of doxorubicin for the treatment of malignant brain tumors in rats. *Anticancer Res* 2005;25:3825–3831.
- [4] 李慧,王春华, 罗 姗,王红冰,黄肖武. 关于小鼠脑内不同注射手法的探讨. *Laboratory Animal Science* 2013 April;Vol.30 No.2.
- [5] AM Sonabend, IV Ulasov, Y Han, CE Rolle, S Nandi, D Cao, MA Tyler, and MS Lesniak. Biodistribution of an oncolytic adenovirus after intracranial injection in permissive animals: a comparative study of Syrian hamsters and cotton rats. *Cancer Gene Ther.* 2009 April; 16(4): 362–372. doi:10.1038/cgt.2008.80.
- [6] Rebecca L. Lowery, Ania K. Intracranial Injection of Adeno-associated Viral Vectors. *Journal of Visualized Experiments.*2010.
- [7] Jian-Ping Zhang, Qi Xu, Xiang-Shan Yuan, Yoan Cherasse, Serge N. Schiffmann, Alban de Kerchove d'Exaerde, Wei-Min Qu, Yoshihiro Urade, Michael Lazarus, Zhi-Li Huang and Rui-Xi Li. Projections of nucleus accumbens adenosine A2A receptor neurons in the mouse brain and their implications in mediating sleep-wake regulation. *Frontiers in Neuroanatomy* www.frontiersin.org December 2013 ,Volume 7,Article 43,1.
- [8] Mi-La Kim, Shengjun Han, Sat-Byol Lee, Jung Hye Kim, Hee Kyung Ahn, Youngbuhm Huh. Evaluation of recombinant adenovirus-mediated gene delivery for expression of tracer genes in catecholaminergic neurons. *Anat Cell Biol*2010 43:157~164.
- [9] Paul A. Lapchak ), Dalia M. Araujo, Dana C. Hilt, Jackie Sheng, Shoushu Jiao. Adenoviral vector-mediated GDNF gene therapy in a rodent lesion model of late stage Parkinson's disease. *Brain Research* 1997 777:153–160.